

# Lovibond® Colour Measurement Tintometer® Group



## Lovibond® Vermelho, Amarelo, Neutro (Fx)

Artigo No: 169630



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Lovibond® Vermelho, Amarelo, Neutro (Fx)

A set of five glass conformance filters for checking Model Fx instruments in Lovibond RYBN scale.

#### Dados Técnicos

Glass Standard Values	Valores aproximados: -
	0.5R 1.3Y 0.1N
	1.3R 3.3Y 0.5N
	2.5R 8.6Y 0.1N
	4.0R 22Y 0.1N
	7.1R 51Y 0.5N

Filters 5

### Texto do Escopo de Entrega

- Lovibond RYBN 0.5R 1.3Y 0.1N filter
- Lovibond RYBN 1.3R 3.3Y 0.5N filter
- Lovibond RYBN 2.5R 8.6Y 0.1N filter
- Lovibond RYBN 4.0R 22Y 0.1N filter
- Lovibond RYBN 7.1R 51Y 0.5N filter
- Certificate of conformity