

Lovibond® Water Testing

Tintometer® Group



Kit de teste, DBNPA, 0-6.8 mg/L

For the standard situation in water testing



- For testing cooling water and other industrial applications
- Test equipment with Checkit® Comparator
- Reagents for up to 100 tests
- Easy to use & result evaluation with visual comparison

Artigo No: 56K701190

A fast biocide

DBNPA (dibromonitrile propionamide) is a fast-acting biocide that rapidly degrades to non-hazardous by-products, including ammonia, bromide ions, dibromoacetonitrile and dibromoacetic acid. It therefore also meets strict environmental requirements. It is one of the biocides recommended by many guidelines for the treatment of industrial process water.

Reliably eliminating microorganisms

In water, DBNPA kills microorganisms quickly and acts reliably against bacteria, fungi, yeasts, algae and microbiological slime: it can quickly penetrate the microbial cell membrane, act on certain protein genes and thus lead to cell death.

Test Kit with comparator

The application of the Test Kit is very easy. After pre-filtration of the water sample, a certain number of reagents is added step by step. Further filtration and the addition of reagent tablets follows. The colour change is used to interpret the results with the help of the colour disc in the comparator.

Indústria

Outras Indústrias

Aplicação

Controle de Desinfecção

Kit de teste, DBNPA, 0-6.8 mg/L

The use of DBNPA in cooling water treatment is widespread because it is a fast and reliable biocide. It is used to combat microorganisms that can cause damage to the water system and material. To ensure that it can work efficiently, regular checks of the DBNPA content should be carried out. This test kit is a reliable companion for these routine tasks.

Faixa de Medição

Test Name	Faixa de Medição
DPNPA	0 - 6.8 mg/L DBNPA

Texto do Escopo de Entrega

- Plastic case
- CHECKIT® Comparator with standard optics
- CHECKIT® Disc DBNPA
- 30 tablet reagents DPD No.1 rapid
- 30 tablet reagents DPD No.3 rapid
- 2 Plastic tube with lid, 10 mL
- Instruction manual

Título	Artigo No
Comparador CHECKIT®	145000
Cubeta de plástico, profundidade de camada de 13,5 mm, com tampa cinzenta, jogo de 10	145500
Cubeta de plástico, profundidade de camada de 13,5 mm, com tampa cinzenta, jogo de 5	145505
Cubeta de plástico, profundidade de camada de 13,5 mm, com tampa cinzenta, jogo de 100	145510
Seringa, plástico, 20 mL	56A006501

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Alemanha

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Reino Unido

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malásia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brasil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
EUA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Índia

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Espanha

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies