

Lovibond[®] Colour Measurement Tintometer[®] Group



Número de Cor do AOCS Tintometer 2.5R 20.0Y (5¼")

Artigo No: 134270



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Número de Cor do AOCS Tintometer 2.5R 20.0Y (5¼")

Garrafa de 500 ml, certificada segundo o sistema QS ISO 9001 com um prazo de validade de 12 meses.

Rastreabilidade completa de acordo com normas reconhecidas internacionalmente

Cada garrafa é fornecida com certificados completos, incluindo ficha de dados de segurança.