

# Lovibond® Colour Measurement Tintometer® Group



## Corante padrão líquido Cobalto platina número de cor 5 - 500 ml

Artigo No: 134140



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Corante padrão líquido Cobalto platina número de cor 5 - 500 ml

Garrafa de 500 ml, certificada segundo o sistema QS ISO 9001 com um prazo de validade de 12 meses.

Rastreabilidade completa de acordo com normas reconhecidas internacionalmente

Cada garrafa é fornecida com certificados completos, incluindo ficha de dados de segurança.