

Lovibond® Colour Measurement Tintometer® Group



Conjunto de filtros EP castanho/amarelo

Artigo No: 139440



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Conjunto de filtros EP castanho/amarelo

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown/Yellow scale.

Dados Técnicos

Glass Standard Values	BY3, BY5
Filters	2

Texto do Escopo de Entrega

- EP Brown/Yellow 3 filter
- EP Brown/Yellow 5 filter
- Certificate of Conformity