

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Lovibond<sup>®</sup> RYBN 色标 2.1R 11.0Y 0.5N (5¼)

货号: 134110



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Lovibond<sup>®</sup> RYBN 色标 2.1R 11.0Y 0.5N (5¼)

500 ml 瓶子，经过 ISO 9001 QS 系统认证，保质期为 12 个月。  
完全可追溯性符合国际公认标准  
每个瓶子均随附完整证书，包括安全数据表。